



MICAP-OES 1000

Direct analysis of trace elements in Middle Distillate Fuels with a compact N₂ ICP-OES system

Introduction

Middle distillates are a category of fuels produced from the distillation of crude oil. This includes products commonly referred to as diesel, marine diesel oil, kerosene and jet fuel. These colorless to light yellow fuels each differs slightly in their components based on the application. Elemental impurities present can impact critical properties such as combustion, corrosion, shelf life, and production of engine deposits. Inductively Coupled Plasma Optical Emission Spectroscopy (ICP-OES) is commonly utilized to monitor the trace elemental content of these various fuels.

This note describes the use of the simultaneous MICAP™ OES-1000 N₂ ICP-OES system for routine analysis of middle distillate fuels following the procedures outlined in ASTM Method D7111. Analytical performance of the MICAP is presented to demonstrate accurate and stable performance during the direct analysis of several fuel samples.

Instrumentation

The ASTM standard test method D7111(1) covers the middle distillate products, which includes distillation fractions between 150 and 390 °C. The samples are analyzed directly against organometallic standards prepared in a kerosene matrix. This requires the Radom Instruments MICAP-OES 1000 to be configured with an organic solvent introduction system. A double-pass cyclonic spray chamber and high-solids V-Groove nebulizer were connected to the standard one-piece torch for introduction of samples into the N₂ plasma.

Microwave energy in the MICAP is coupled into the Cerawave™ ring in a highly efficient process that creates the magnetic fields required to inductively couple energy into the robust N₂ plasma source. The 4 MP sCMOS camera simultaneously collects the resulting emission lines from the high-resolution spectrometer.

Experimental Conditions

The analyses were performed on the simultaneous MICAP (Figure 1) utilizing the conditions listed in Table 1. Samples were introduced with an autosampler, with each sample analysis taking 3 minutes to complete.

The MICAP in this work was uniquely configured with air (rather than the standard N₂) as the nebulizer gas. This modification reduces the continuum background created by the combustion of carbon, resulting in significantly improved detection limits for several elements. Table 2 provides specifics on the analyte and Yttrium (Y) internal standard (I/S) wavelengths used for this analysis, as well as the detection limits obtained from 3x the standard deviation of 10 blank measurements.



Figure 1. MICAP-OES 1000 N₂ ICP-OES

Table 1. Operational conditions

Parameter	Value
Torch	Quartz 1-piece, 1.5mm injector
Spray Chamber	Double pass cyclonic
Nebulizer	V-Groove glass
Sample Tubing	Solvaflex Blk/Blk (0.76 mm ID)
Drain Tubing	Solvaflex Blu/Yel (1.52 mm ID)
Rinse	PremiSolv™
Stabilization Time	15 sec
Power	1000 W
Coolant Gas Flow	16 L/min N ₂
Auxiliary Gas Flow	0.4 L/min N ₂
Nebulizer Gas Flow	0.55 L/min Air
Peristaltic Pump	90 rpm
Plasma Viewing	Axial
Camera Exposure	20 sec (125 ms @ 160 exps)
# of Repeats	3

Table 2. Analyte wavelengths and fuel detection limits

Analyte Wavelength (nm)	Detection Limit (mg/Kg, ppm)
Ag 328.068	0.001
Al 396.152	0.008
B 249.677	0.032
Ba 585.368	0.024
Ca 393.366	0.001
Cd 228.802	0.008
Cr 428.973	0.004
Cu 324.754	0.0008
Fe 259.940	0.014
K 766.490	0.042
Mg 279.553	0.0007
Mn 257.610	0.003
Mo 281.615	0.021
Na 588.995	0.008
Ni 300.249	0.015
P 253.560	0.170
Pb 283.305	0.028
Si 251.611	0.006
Sn 283.998	0.015
Ti 334.940	0.0008
V 309.310	0.008
Zn 213.857	0.023
Y 371.029 (I/S)	-

Standard and Sample Preparation

The calibration standards, blanks, and samples were prepared in accordance with the protocols outlined in ASTM Method D7111. Working standards were prepared by diluting oil-based stock standards (VHG, LGC Standards, Manchester, NH, USA) in PremiSolv™ (Conostan, AnalytiChem, Bale-D'Urfe, QC, Canada) prepared by weight. The Yttrium (Y) internal standard was added at a final concentration of 3 mg/Kg (ppm). A blank and a 2 mg/Kg mixed standard were prepared for all the listed analytes. Additionally, a 5 mg/Kg Phosphorus (P) standard was also prepared to provide a higher P calibration range.

Three middle distillate samples were sourced locally. These included a diesel fuel, an off-road (agricultural) diesel fuel, and Jet A fuel sample. Each of the samples were also prepared with a known spike (at 1 mg/Kg) to evaluate system performance. The Y internal standard was also added to each of these samples.

Results

The MICAP was calibrated as described above and all wavelengths were viewed in the Radom Instruments Software (RIS) Profiles View to verify no spectral interferences existed and to set peak integration points. The emission lines utilized for quantitation were selected from numerous line options available. The simultaneous, full wavelength range data collection also provided the ability during method development to confirm accurate measurements across multiple emission lines for each element.

Figure 2 provides example emission peaks for the analytes of Mg, Cu, and Ti. Note the clear baseline resolution between the Ti 334.903 and 334.940 nm emission lines, demonstrating the high resolving power of the spectrometer. This capability helps ensure the analyte peaks of interest are not impacted by nearby spectral interferences from other potential elemental components.

The middle distillate fuel results and analytical spike (at 1 mg/Kg) results are displayed in Table 3. The three samples show very little elemental contamination is present. The sample spike results show excellent recoveries, with the exception of Boron (B) - for which stability / volatility issues are quite common.

Internal standard recovery values were consistently within $\pm 15\%$, indicating that the system response was very stable between the calibration standards and the samples.




Figure 2. Profiles view of Mg, Cu & Ti emission lines in 1 mg/Kg spiked Jet A fuel. Sample emission lines displayed in brown with blank baseline indicated in green.

Table 3. Middle Distillate Fuel Analysis - Sample and 1 mg/Kg Spike Results

Analyte Wavelength (nm)	Diesel Fuel		Agricultural Diesel Fuel		Jet A Fuel	
	Conc. (mg/Kg)	Spike Recovery (%)	Conc. (mg/Kg)	Spike Recovery (%)	Conc. (mg/Kg)	Spike Recovery (%)
Ag 328.068	-0.002	97	-0.001	97	0.000	91
Al 396.152	0.000	98	-0.002	98	-0.021	92
B 249.772*	0.100	86	0.072	93	0.062	63
Ba 585.368	0.088	98	0.085	99	-0.239	98
Ca 393.366	-0.003	108	0.015	107	0.009	102
Cd 228.802	0.019	100	0.013	99	-0.021	95
Cr 428.973	-0.012	97	-0.013	96	0.003	91
Cu 324.754	0.000	100	0.000	99	0.001	93
Fe 259.940	0.030	98	0.094	100	-0.013	98
K 766.490	0.101	90	0.114	88	-0.064	95
Mg 279.553	0.000	99	0.009	98	0.002	98
Mn 257.610	0.007	97	0.010	97	-0.011	96
Mo 281.615	-0.013	96	-0.013	98	0.005	96
Na 588.995	0.035	92	0.039	92	-0.078	90
Ni 300.249	-0.044	98	-0.001	99	-0.000	92
P 253.560	-0.342	94	-0.175	100	0.821	95
Pb 283.305	0.009	100	0.022	97	-0.051	90
Si 251.611	0.007	99	0.042	99	0.044	92
Sn 283.998	-0.025	97	-0.002	98	-0.023	91
Ti 334.940	-0.006	98	-0.004	98	0.008	98
V 309.310	-0.011	97	-0.012	97	0.016	96
Zn 213.857	-0.022	100	-0.022	101	-0.045	94

*Boron spike results showed some variability, likely due to common stability/volatility issues with this element

Conclusions

The MICAP-OES 1000 N₂ ICP-OES demonstrated its capability for the direct analysis of middle distillate fuels in accordance with ASTM Standard Method D7111. No carbon deposition from the samples was observed on the MICAP torch throughout this application work, demonstrating a significant impact on productivity as no user maintenance is required to clean the torch and then restart the system. The excellent accuracy and this stability confirms MICAP's capability for routine analysis of these fuel samples.

References

1. ASTM D7111, Standard Test Method for Determination of Trace Elements in Middle Distillate Fuels by Inductively Coupled Plasma Atomic Emission Spectrometry (ICP-AES), ASTM Intn, West Conshohocken, PA, 2025, www.astm.org
2. ASTM D5708, Standard Test Method for Determination of Nickel, Vanadium, and Iron in Crude Oils and Residual Fuels by Inductively Coupled Plasma (ICP) Atomic Emission Spectrometry, ASTM Intn, West Conshohocken, PA, 2025, www.astm.org

Scan the code to access our website for additional information and our resources

